

Product Description

AcuGenix™ 1st cDNA Synthesis Premix for qPCR (+gDNA wiper) contains all reagents required for first-strand cDNA synthesis and is suitable for two-step qRT-PCR detection. 4×gDNA Wiper Mix can thoroughly eliminate residual genomic DNA contamination in RNA templates, ensuring more accurate quantification results. It also simplifies qPCR primer design, eliminating the need for intron-spanning primers. The 5×Neoscript 1st cDNA Synthesis Premix is an efficient and convenient first-strand cDNA synthesis master mix. The reverse transcription reaction can be performed simply by adding the RNA template and nuclease-free water.

cDNA synthesized with this product is compatible with both fluorescent dye and probe-based qPCR platforms. Corresponding supporting reagents can be selected according to experimental purposes for high-performance gene expression analysis.

Components

Components	BR1D301-01 100 T (20 µL/T)	BR1D301-02 500 T (20 µL/T)
4×gDNA Wiper Mix	40 µL	400 µL
5×Neoscript 1st cDNA Synthesis Premix ¹	40 µL	400 µL
5×No RT Control Mix ²	4 µL	40 µL
RNase-free ddH ₂ O	200 µL	2×1mL

1. This reagent contains components that inhibit DNase activity. Samples treated with 4×gDNA Wiper Mix can be directly subjected to reverse transcription for cDNA synthesis, enabling rapid completion of the entire process from genomic DNA removal to cDNA synthesis. It contains Random 6 mers, Oligo (dT)18 Primer and dNTPs, allowing efficient and uniform synthesis of various cDNAs.

2. Excluding Neoscript RTase, this reagent shares the same composition as 5×Neoscript 1st cDNA Synthesis Premix, and is used for preparing No RT control reactions.

Storage

Store at -20±5°C.

Notes

- For Research Use Only. Not for use in diagnostic procedures.
- Mix thoroughly before use and avoid repeated freeze-thaw cycles.
- For long-fragment amplification, to prevent destruction of the cDNA structure, incubation at 70 °C for 15 min is recommended for enzyme inactivation.
- The addition of RNase Inhibitor can effectively inhibit RNA template degradation caused by nucleases introduced during RNA extraction.
- High-purity and intact RNA templates are critical for reverse transcription to generate high-quality and full-length cDNA strands. The addition of gDNA Wiper maximizes the accuracy and reproducibility of qPCR results.

Prepare RT Reaction Mix

1. Reaction system setup (prepare on ice).

Components	Volume per Reaction
4×gDNA Wiper Mix	5 µL
Template RNA	10 pg -1 µg*
RNase-free ddH ₂ O	To 16 µL

*Add an appropriate amount of RNA template according to experimental requirements.

2. Mix thoroughly by pipetting up and down, then incubate at 42°C for 2 minutes (temperature control can be performed using a thermal cycler).

3. Add the following components to the above system to a final volume of 20 µL (prepare on ice), and mix gently.

Components	Volume per Reaction
5×Neoscript 1st cDNA Synthesis Premix	4 µL

4.Reverse transcription program: 50 °C for 15 min; 85 °C for 5 s.

*For templates with complex secondary structures or high-GC regions, the reaction temperature can be appropriately increased.

5.The reverse transcription product can be directly used for conventional PCR and real-time fluorescent quantitative PCR, or stored at -20 °C.

Prepare Reaction Mix

The following is the operational procedure for subsequent qPCR detection on a real-time fluorescent quantitative PCR instrument using AcuGenix™ FastAmpli qPCR Master Mix (Universal ROX) (Cat.: BR1E102). If used with other Biori qPCR products, please refer to the corresponding product manual.

Components	Volume per Reaction
2×AcuGenix™ FastAmpli qPCR Master Mix (Universal ROX)	10 µL
Forward Primer (10 µM)	0.4 µL
Reverse Primer (10 µM)	0.4 µL
cDNA*	-
RNase-free ddH ₂ O	To 20 µL

*It is recommended to dilute cDNA with RNase-free water, and the volume of added cDNA shall not exceed 1/10 of the total qPCR reaction system volume.

Reaction Program

Conventional PCR procedures			
Steps	Temp	Time	Cycles
Initial denaturation ¹	95 °C	1 min	1
Degeneration ²	95 °C	10 s	40
*Annealing and Elongation ³	56-64 °C	30 s	
*Melting curve analysis ⁴			

Fast PCR procedure			
Steps	Temp	Time	Cycles
Initial denaturation ¹	95 °C	30 s	1
Degeneration ²	95 °C	1 s	40
*Annealing and Elongation ³	56-64 °C	10 s	
*Melting curve analysis ⁴			

1.The pre-denaturation conditions are applicable to most amplification reactions. The conventional PCR procedure uses 1 min, while the fast PCR procedure can be set to a minimum of 30 s. For templates with complex structures, the pre-denaturation time can be extended to 3 minutes to improve the pre-denaturation effect.

2.Denaturation: 10 s for the conventional PCR procedure, and a minimum of 1 s for the fast PCR procedure.

3.Annealing/Extension: 30 s for the conventional PCR procedure. For the fast PCR procedure, the extension time can be set to 10 s for amplicons within 200 bp. For products longer than 200 bp, it is recommended to appropriately extend the extension time or adopt a three-step protocol. The annealing and extension temperature shall be adjusted according to the T_m value of the designed primers.

4.Melting curve analysis: Please set the program recommended by the fluorescence quantitative PCR instrument used.

5.* Set signal acquisition at marked position.