

Product Description

T7 RNA Polymerase GMP-grade is an in-house developed and multi-step purified recombinant T7 RNA polymerase manufactured by Biori Biotech to GMP standards. It exhibits high transcriptional catalytic activity across diverse template types and nucleotide substrates. Under appropriate reaction buffer conditions, this enzyme enables efficient synthesis of large quantities of RNA from DNA templates. T7 RNA Polymerase GMP-grade specifically recognizes the T7 promoter sequence (5'-TAATACGACTCACTATAG-3') and initiates transcription from the G residue at this site, converting downstream DNA sequences into single-stranded RNA. Using natural or modified nucleotides as substrates, a single 20 μ L in vitro transcription reaction can yield more than 200 μ g of RNA.

Components

Components	Cat. No.	Quantity	Volume
T7 RNA Polymerase GMP-grade (50 U/ μ L)	GMP-BP-E01-5K	5 KU	100 μ L
	GMP-BP-E01-50K	50 KU	1 mL
	GMP-BP-E01-500K	500 KU	10 mL

Storage

Store at $-20\pm 5^{\circ}\text{C}$.

Product Information

Product Name	T7 RNA Polymerase GMP-grade
Source	Recombinant <i>E.coli</i>
Activity	50 U/ μ L
Unit Definition	One unit (U) is defined as the amount of enzyme required to incorporate 1 nmol of NTP into RNA in 1 hour at 37°C , pH 8.0.
Storage Buffer	50 mM Tris-HCl, 1 mM EDTA, 10 mM DTT, 100 mM NaCl, 0.1%(v/v) Triton X-100, 50%(v/v) glycerol, pH7.9 at 25°C

Quality control

1. Solution appearance: clear and transparent, free of visible particulate matter.
2. Activity >50 U/ μ L.
3. Protein purity $\geq 95\%$.
4. Free of exogenous DNase, RNase, exonuclease, and endonuclease activity.
5. Residual host-cell DNA: ≤ 100 pg/mg.
6. Residual host-cell protein: <50 ppm.
7. Heavy metals <10 ppm.
8. HBV, HCV, HIV, and mycoplasma: not detected.
9. Bacterial endotoxin: <5 EU/mL.
10. pH 7.0-8.0.

Recommended Transcription System

Components	Volume
RNase free ddH ₂ O	To 20 μ L
ATP Solution(100 mM)	2 μ L
CTP Solution(100 mM)	2 μ L
GTP Solution(100 mM)	2 μ L
UTP Solution(100 mM)	2 μ L
10 \times Transcription Buffer	2 μ L
RNase Inhibitor GMP-grade (40 U/ μ L)	1 μ L
Pyrophosphatase, Inorganic GMP-grade (0.1 U/ μ L)	1 μ L
T7 RNA Polymerase GMP-grade(50 U/ μ L)	2 μ L
Template DNA	1 μ g

*Incubate at 37°C for 2 hours.

Notes

1. In vitro transcription reactions are highly sensitive to RNase. Strictly avoid introducing RNase into the reaction system. All laboratory consumables, including pipette tips and microcentrifuge tubes, must be RNase-free.
2. Equilibrate the 10×Transcription Buffer to room temperature before use. At low temperatures, DNA may co-precipitate with spermidine, reducing transcription yield; the DNA template should therefore always be added last. Recommended addition order: nuclease-free water → NTPs → 10×Transcription Buffer → RNase Inhibitor → Pyrophosphatase, Inorganic → T7 RNA Polymerase → Template DNA.
3. Wear gloves throughout the experimental procedure.
4. Mix thoroughly before use. Avoid repeated freeze-thaw cycles.