

## Product Description

DNase I GMP-grade(Deoxyribonuclease I) is an endonuclease that efficiently digests single-stranded and double-stranded DNA. DNase I requires  $\text{Ca}^{2+}$  for activity and is activated by  $\text{Mg}^{2+}$  or  $\text{Mn}^{2+}$ . In the presence of  $\text{Mg}^{2+}$ , DNase I cleaves DNA at random sites; in the presence of  $\text{Mn}^{2+}$ , the enzyme cleaves both strands of double-stranded DNA at approximately the same position, generating blunt-ended fragments or fragments with 1-2 nucleotide overhangs. DNase I exhibits optimal endonuclease activity at pH 7-8 and is suitable for removing single-stranded or double-stranded DNA from RNA preparations.

## Components

Components	GMP-BP-E04-200	GMP-BP-E04-2K	GMP-BP-E04-20K
	200 U	2 KU	20 KU
DNase I GMP-grade (2 U/ $\mu\text{L}$ )	0.1 mL	1 mL	10 mL
10 $\times$ DNase I Reaction Buffer	0.2 mL	2 mL	20 mL

## Storage

Store at  $-20\pm 5^\circ\text{C}$ .

## Product Information

Product Name	DNase I GMP-grade
Source	Recombinantly expressed in Pichia yeast strain
Activity	2 U/ $\mu\text{L}$
Unit Definition	One unit (U) is defined as the amount of enzyme required to completely digest 1 $\mu\text{g}$ of plasmid DNA in 10 min at $37^\circ\text{C}$ .
Storage Buffer	10 mM Tris-HCl, 2 mM $\text{CaCl}_2$ , 50% (v/v) glycerol, pH 7.6.
Inactivation	Add EDTA to a final concentration of 5 mM and incubate at $65^\circ\text{C}$ for 10 min to inactivate the enzyme.

## Quality control

1. Solution appearance: clear and transparent, free of visible particulate matter.
2. Activity >2 U/ $\mu\text{L}$ .
3. Protein purity  $\geq 95\%$ .
4. Free of exogenous RNase activity.
5. Residual host-cell DNA:  $\leq 100$  pg/mg.
6. Residual host-cell protein: <50 ppm.
7. Heavy metals <10 ppm.
8. HBV, HCV, HIV, and mycoplasma: not detected.
9. Bacterial endotoxin: <5 EU/mL.
10. pH 7.8-8.0

## Recommended Usage (DNA Removal from RNA Samples)

Using RNase-free microcentrifuge tubes or PCR tubes, prepare the reaction mixture as follows:

Components	Volume
10 $\times$ DNase I Reaction Buffer	1 $\mu\text{L}$
DNase I GMP-grade (2 U/ $\mu\text{L}$ )	1 $\mu\text{L}$
RNA	Up to 1 $\mu\text{g}$
RNase-free ddH <sub>2</sub> O	To 10 $\mu\text{L}$

Mix gently and incubate at  $37^\circ\text{C}$  for approximately 30 min.

## Application

1. Preparation of DNA-free RNA samples for downstream RNA analysis.
2. Digestion of linearized plasmid template following in vitro transcription.
3. Removal of residual genomic DNA prior to RT-PCR.
4. DNase I footprinting analysis of DNA-protein interactions.

**Notes**

1. Mix the reaction mixture by gentle pipetting; avoid vigorous vortexing.
2. Keep the enzyme on ice throughout the procedure.
3. Mix thoroughly before use. Aliquot for storage to avoid repeated freeze-thaw cycles.
4. The optimal amount of DNase I GMP grade should be optimized according to the specific experimental requirements.