

## Product Description

AcuGenix™ Taq DNA Polymerase is a hot-start Taq that not only better suppresses non-specific reactions caused by primer false annealing or primer dimers during PCR system preparation and amplification, making the product highly specific and more effective for amplifying low-concentration templates, suitable for multiplex PCR amplification reactions; but also has good system adaptability, achieving stable amplification results in different types of PCR reactions.

## Components

Components	BR1A101-01 100 U	BR1A101-02 1000 U	BR1A101-03 5000 U
5 U/μL AcuGenix™ Taq DNA Polymerase	0.02 mL	0.2 mL	1 mL

## Unit Definition

One unit of activity (U) refers to the amount of enzyme needed to incorporate 10 nmol of deoxyribonucleotides into acid-insoluble material using activated salmon sperm DNA as a template/primer within 30 minutes at 74°C.

## Storage

Storage at -20±5°C.

## Notes

- 1.For Research Use Only. Not for use in diagnostic procedures.
- 2.Mix thoroughly before use and avoid repeated freeze-thaw cycles.
- 3.This product can be quickly hot-started at 95°C or 94°C for 1-5 minutes.
- 4.This system has strong adaptability and higher specificity and sensitivity.
- 5.It is suitable for high-sensitivity PCR detection reagents and can be used in multiplex PCR amplification reactions.
- 6.This product has both 5'-3' polymerase and exonuclease activity, no 3'-5' exonuclease activity or proofreading function.
- 7.This product is suitable for qualitative/quantitative detection in regular PCR and RT-PCR.
- 8.PCR products have an A-tail at the end, which can be directly cloned into T vector.
- 9.For primers with low annealing temperature or amplification of fragments longer than 200 bp, it is recommended to use a three-step method.

## Quality Control

- 1.SDS-PAGE electrophoresis purity no less than 98%.
- 2.Amplification sensitivity, batch-batch difference, and stability.
- 3.No exogenous nuclease activity, no exogenous endonuclease or exonuclease contamination.

## Prepare Reaction Mix

Component	Volume per Reaction	Final Concentration
10×PCR Buffer II (Mg <sup>2+</sup> free) <sup>1</sup>	5 μL	1×
dNTPs (10 mM each)	1 μL	200 μM
25 mM MgCl <sub>2</sub>	2-8 μL	1-4 mM
5 U/μL AcuGenix™ Taq DNA Polymerase	0.25-0.5 μL	1.25-2.5U
25×Primer Mix <sup>2</sup>	2 μL	1×
Template	--	<1 μg/Reaction
ddH <sub>2</sub> O	To 50 μL	--

- 1.This buffer does not contain dNTPs or Mg<sup>2+</sup>, so they must be added to the reaction system before use.
- 2.If used for qPCR/qRT-PCR, a fluorescent probe needs to be added to the reaction system. Typically, a final primer concentration of 0.2 μM performs good results. If the reaction performance is poor, adjust the primer concentration within the range of 0.2-1 μM. Typically, probe concentrations are optimized in the range of 0.1-0.3 μM. Combinations of primers and probes can be tested using gradient experiments to find their optimal combination.

## Reaction Program

Two-step Method			
Procedure	Temp	Time	Cycles
Initial denaturation	95°C	1-5 min	1
Degeneration	95°C	10-20 s	35-50
Annealing and Elongation	56-64°C	20-60 s	
-			

Three-step Method			
Procedure	Temp	Time	Cycles
Initial denaturation	95°C	1-5 min	1
Degeneration	95°C	10-20 s	35-50
Annealing	56-64°C	10-30 s	
Elongation	72°C	10-60 s	