

Product Description

AcuGenix™ Reverse Transcriptase is a reverse transcriptase derived from the M-MLV gene of Moloney Murine Leukemia Virus via mutation screening, and expressed in *Escherichia coli*. This enzyme lacks RNase H activity and exhibits higher temperature tolerance. It is suitable for high-temperature reverse transcription, which helps eliminate the adverse effects of RNA advanced structures and non-specific factors on cDNA synthesis, and possesses higher stability and reverse transcription synthesis capacity.

Components

Components	BR1D101-01 10,000 U	BR1D101-02 40,000 U	BR1D101-03 200,000 U
200 U/μL AcuGenix™ Reverse Transcriptase	0.05 mL	0.2 mL	1 mL

Unit Definition

One unit of activity (U) refers to the amount of enzymes required to incorporate 1 nmol dTTP into poly(A)•Oligo(dT)25 template/primer at 37°C within 10 minutes.

Storage

Storage at -20±5°C.

Notes

- 1.This product is suitable for temperature optimization of reverse transcription between 42°C-55°C.
- 2.This product has better stability, suitable for high-temperature reverse transcription amplification, effectively passing through RNA complex structure regions, applicable to one-step multiple fluorescent quantitative RT-PCR detection.
- 3.This product is compatible with various PCR amplification enzymes and suitable for highly sensitive RT-PCR reactions.
- 4.This product is suitable for highly sensitive one-step fluorescent quantitative RT-PCR reaction, effectively improving the detection rate of low-concentration templates.
- 5.This product is suitable for cDNA library construction.
- 6.This product is suitable for 3' and 5' RACE.

Quality Control

- 1.SDS-PAGE electrophoresis purity no less than 98%.
- 2.Amplification sensitivity, batch-batch difference, and stability.
- 3.No exogenous nuclease activity, no exogenous endonuclease or exonuclease contamination.

First-strand synthesis

- 1.Prepare the reaction system according to the table below;

Component	Quantity per Reaction
Oligo (dT) ₁₂₋₁₈ Primer	50 pmol
Or Random Primer*	50 pmol 20-100 pmol
Or Gene specific Primer*	2 pmol
dNTPs (10 mM each)	1 μL
Template RNA	total RNA≤5 μg, mRNA≤1 μg
RNase-free ddH ₂ O	To 10 μL

*Select different types of primers according to experimental needs.

- 2.Heat at 65°C for 5 minutes, then quickly cool on ice for 2 minutes.
- 3.Add the following components to the above system until the total volume reaches 20 μL, and gently mix.

Component	Volume per Reaction
5×First-Strand Buffer	4 μL
200 U/μL AcuGenix™ Reverse Transcriptase	1 μL
40 U/μL AcuGenix™ RNase Inhibitor	1 μL
RNase-free ddH ₂ O	To 20 μL

4. Perform reactions under the following conditions:

(1) If using Random Primer random primers, incubate at 25°C for 10 minutes and then at 50°C for 30-60 minutes.

- (2) If using Oligo dT or specific primers, incubate at 50°C for 30-60 minutes.
5. Heat inactivate AcuGenix™ Reverse Transcriptase by incubating at 95°C for 5 minutes and terminate the reaction.
6. The reverse transcription product can be directly used in PCR reactions and fluorescence quantitative PCR reactions or stored long-term at -20°C.

Prepare Reaction Mix

Components	Volume per Reaction	Final Concentration
10×PCR Buffer II (Mg ²⁺ free) ¹	5 µL	1×
dNTPs (10 mM each)	1 µL	200 µM
25 mM MgCl ₂	2-8 µL	1-4 mM
5 U/µL AcuGenix™ Taq DNA Polymerase	0.25-0.5 µL	1.25-2.5U/Reaction
25×Primer Mix ²	2 µL	1×
Template	—	1 µg/Reaction
ddH ₂ O	To 50 µL	—

1. 10×PCR Buffer II (Mg²⁺ free) does not contain dNTP or Mg²⁺. Please add dNTPs and MgCl₂ when preparing reaction system.

2. If used for qPCR/qRT-PCR, fluorescent probes need to be added to the reaction system. The final concentration of primers is usually 0.2 µM for better results. When the reaction performance is poor, the primer concentration can be adjusted within the range of 0.2-1 µM. The probe concentration is usually optimized in the range of 0.1-0.3 µM. Concentration gradient experiments can be performed to find the best combination of primers and probes.

Reaction Program

Standard PCR Program			
Step	Temperature	Duration	Cycles
Initial denaturation	95°C	1-5 min	1
Degeneration	95°C	10-20 s	30-40
Annealing and Elongation	50-60°C	10-30 s	
Extension	72°C	10-60 s	