

Product Description

Robustart Tth DNA Polymerase is a hot-start Tth enzyme that has 5'-3' DNA template-dependent polymerase activity, making it suitable for various PCR reactions. It also has higher reverse transcriptase activity in the presence of Mn^{2+} , effectively eliminating the adverse effects of RNA secondary structures and non-specific primer annealing at low temperatures during cDNA synthesis, thereby increasing reverse transcription specificity and improving amplification efficiency and yield. It can be used for single-enzyme one-step RT-PCR reactions. Robustart Tth's activity is blocked at room temperature, which improves its specificity and significantly enhances the amplification effect of low-concentration templates.

Components

1. 5 U/ μ L Robustart Tth DNA Polymerase
2. 10 $\times$ Tth RT Buffer Basic (Mn^{2+} free) (optional)
3. 25 mM Mn (OAc)₂ (optional)

*10 $\times$ Tth RT Buffer Basic (Mn^{2+} free) does not contain dNTPs or Mg^{2+} . Please add dNTPs and Mn(OAc)₂ when preparing reaction systems.

Unit Definition

One unit of activity (U) refers to the amount of enzyme required to incorporate 10 nmol of deoxyribonucleotides into acid-insoluble material using activated salmon sperm DNA as a template/primer within 30 minutes at 74°C.

Storage

Store at $-20\pm 5^{\circ}\text{C}$. Vortex before use. Avoid freeze-thaw cycles.

Quality Control

1. SDS-PAGE electrophoresis purity no less than 98%.
2. Amplification sensitivity, between-batch difference, and stability.
3. No exogenous nuclease activity, no exogenous endo/exonuclease contamination.

Prepare Reaction Mix

Components	Volume per Reaction	Final Concentration
10 $\times$ Tth RT Buffer Basic (Mn^{2+} free) ¹	5 μ L	1 $\times$
dNTPs (10 mM each)	1.5 μ L	300 μ M
25 mM Mn(OAc) ₂	2-6 μ L	1-3 mM
5 U/ μ L Robustart Tth DNA Polymerase	1-2 μ L	5-10U/Reaction
25 $\times$ Primer Mix ²	2 μ L	1 $\times$
Template	--	<1 μ g/Reaction
ddH ₂ O	To 50 μ L	--

1. This Buffer does not contain dNTP or Mn^{2+} . dNTPs and Mn(OAc)₂ must be added before use.
2. If used for qPCR/qRT-PCR, fluorescent probes need to be added to the reaction system. The final concentration of primers is usually 0.2 μ M for better results. when the reaction performance is poor, the primer concentration can be adjusted within the range of 0.2-1 μ M. The probe concentration is usually optimized in the range of 0.1-0.3 μ M. Concentration gradient experiments can be performed to find the best combination of primers and probes.

Reaction Programe

Standard PCR Program			
Procedure	Temp	Time	Cycle
Hot Start	95°C	1-5 min	1
Denaturation	95°C	10-20 s	35-50
Annealing and Extension	56-64°C	20-60 s	

One-step RT-PCR Program			
Procedure	Temp.	Time	Cycle
Denaturation	90°C	30 s	1
Reverse Transcription	60°C	25 min	1
Denaturation	95°C	1-5 min	1
Denaturation	95°C	10-20 s	35-50
Annealing and Extension	56-64°C	20-60 s	

Notes

1. Robustart Tth DNA Polymerase has 5'-3' polymerase and exonuclease activities but no 3'-5' exonuclease activity. PCR products end with A.
2. The commonly used reverse transcription temperature for Robustart Tth is 60°C, which can be optimized between 60°C -70°C depending on the characteristics of the amplification reaction. Reverse transcription time can be optimized between 15-30 min.
3. It is more suitable for RT reaction with specific primers, and the T_m value of the primer should be 60°C or higher. Primers with low T_m value such as Oligo (dT)₁₈₋₂₀ or Random Primers are not recommended.
4. It can effectively reduce the low reverse transcription efficiency caused by the complex secondary structure of RNA, and improve the specificity of primer-template hybridization.
5. Compared with M-MLV, it has higher thermal stability. And it has higher tolerance to various PCR inhibitors.
6. Due to the use of Mn²⁺, the fidelity of amplification using Tth is reduced. So it is not suitable for cloning, sequencing and other experiments with high fidelity requirements.
7. Three-step PCR method is recommended for primers with low annealing temperature or for the amplification of long fragments over 200 bp.