

## Product Description

10× Tris-Glycine-SDS Running Buffer is commonly used as the electrophoresis buffer for protein SDS-PAGE. This product is supplied as a 10× concentrate and only needs to be diluted with deionized water to a 1× working solution before use.

## Components

Components	BR4C361-01
10× Tris-Glycine-SDS Running Buffer	500 mL

## Storage

Store at 10-30°C.

## Notes

1. For Research Use Only. Not for use in diagnostic procedures.
2. If crystals form, this is normal. Dissolve completely in a 40°C water bath and mix thoroughly before use.
3. For your safety and health, please wear a lab coat and disposable gloves when operating.